

Anti-Rat BCL2L15 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02016

3th Edition

Description:BCL2L15 (BCL2 Like 15) is a Protein Coding gene. Among its related pathways are Regulation of Apoptosis by Parathyroid Hormone-related Protein.

Antigen:Recombinant human Bcl-2-like protein 15 protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: