

Anti-Mouse PDGFRA Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02023

3th Edition

Description: Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

Antigen: Recombinant human PDGFRA protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in primary and metastatic colon tumors and in normal colon tissue. Tumors may express a different isoform to that found in normal tissue.