

**Anti-Mouse PDGFRA Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02023**

---

3th Edition

**Description:** Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

**Antigen:** Recombinant human PDGFRA protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in primary and metastatic colon tumors and in normal colon tissue. Tumors may express a different isoform to that found in normal tissue.