

**Anti-Mouse/Rat SERPINE1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02026**

---

3th Edition

**Description:** Serine protease inhibitor. This inhibitor acts as 'bait' for tissue plasminogen activator, urokinase, protein C and matriptase-3/TMPRSS7. Its rapid interaction with PLAT may function as a major control point in the regulation of fibrinolysis.

**Antigen:** Recombinant human SERPINE1 protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Found in plasma and platelets and in endothelial, hepatoma and fibrosarcoma cells.