

Anti-Human/Mouse SEPT5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02028

3th Edition

Description: Filament-forming cytoskeletal GTPase (By similarity). May play a role in cytokinesis (Potential). May play a role in platelet secretion (By similarity).

Antigen: Recombinant human SEPT5 protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed at high levels in the CNS, as well as in heart and platelets (at protein level).