

**Anti-Human/Mouse ANXA10 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02037**

---

3th Edition

**Description:** This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined.

**Antigen:** Synthesized peptide derived from the C-terminal region of human Annexin X

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**