

Anti-Human ANT3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02038

3th Edition

Description: This gene is a member of the mitochondrial carrier subfamily of solute carrier protein genes. The product of this gene functions as a gated pore that translocates ADP from the cytoplasm into the mitochondrial matrix and ATP from the mitochondrial matrix into the cytoplasm. The protein is implicated in the function of the permeability transition pore complex (PTPC), which regulates the release of mitochondrial products that induce apoptosis. The human genome contains several non-transcribed pseudogenes of this gene. SLC25A6 (Solute Carrier Family 25 Member 6) is a Protein Coding gene. Diseases associated with SLC25A6 include Influenza and Autosomal Dominant Progressive External Ophthalmoplegia. Among its related pathways are Respiratory electron transport, ATP synthesis by chemiosmotic coupling, and heat production by uncoupling proteins. and Integration of energy metabolism. GO annotations related to this gene include transporter activity and ATP:ADP antiporter activity. An important paralog of this gene is SLC25A5.

Antigen: Synthesized peptide derived from the Internal region of human ANT3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: