

Anti-Human/Mouse/Rat ALPL Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02061

3th Edition

Description: A phosphate monoester + H₂O = an alcohol + phosphate.

Antigen: Recombinant human ALpp protein

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: