

Anti-Human AQP5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02075

3th Edition

Description: Aquaporin-5 is a protein that in humans is encoded by the AQP5 gene. Aquaporin 5 (AQP5) is a water channel protein. Aquaporins are a family of small integral membrane proteins related to the major intrinsic protein (MIP or AQP0). Aquaporin 5 plays a role in the generation of saliva, tears and pulmonary secretions. AQP0, AQP2, AQP5, and AQP6 are closely related and all map to 12q13.

Antigen: Synthesized peptide derived from the C-terminal region of human AQP5.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: