

Anti-Human/Mouse/Rat NCOA4 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02078

3th Edition

Description:Enhances the androgen receptor transcriptional activity in prostate cancer cells. Ligand-independent coactivator of the peroxisome proliferator-activated receptor (PPAR) gamma.

Antigen:Synthesized peptide derived from the Internal region of human ARA70

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Widely expressed. Also detected in adipose tissues and in different cell lines. Isoform Beta is only expressed in testis.