

Anti-Human SLC25A12 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02081

3th Edition

Description: Calcium-dependent mitochondrial aspartate and glutamate carrier. May have a function in the urea cycle.

Antigen: Synthesized peptide derived from the Internal region of human ARALAR

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: High levels in heart and skeletal muscle, low in brain and very low in kidney.