

Anti-Human/Mouse/Rat ARMC6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02091

3th Edition

Description: The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function.

Antigen: Synthesized peptide derived from the C-terminal region of human ARMC6

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: