

**Anti-Human/Mouse/Rat ARSG Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02103**

---

3th Edition

**Description:** Displays arylsulfatase activity at acidic pH with pseudosubstrates, such as p-nitrocatechol sulfate and also, but with lower activity, p-nitrophenyl sulfate and 4-methylumbelliferyl sulfate.

**Antigen:** Synthesized peptide derived from the Internal region of human Arylsulfatase G

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed, with very low expression in brain, lung, heart and skeletal muscle.