



Anti-Mouse/Rat Cleaved-CASP6 p18 (D162) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02105

3th Edition

Description:Involved in the activation cascade of caspases responsible for apoptosis execution. Cleaves poly(ADP-ribose) polymerase in vitro, as well as lamins. Overexpression promotes programmed cell death.

Antigen:Synthesized peptide derived from the Internal region of mouse Caspase-6 p18

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: