

Anti-Human Cleaved-C1R LC (I464) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02115

3th Edition

Description: C1r B chain is a serine protease that combines with C1q and C1s to form C1, the first component of the classical pathway of the complement system.

Antigen: Synthesized peptide derived from the Internal region of human C1r LC

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: