

Anti-Human/Monkey Cleaved-Cathepsin D LC (G65) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02128

3th Edition

Description: Acid protease active in intracellular protein breakdown. Involved in the pathogenesis of several diseases such as breast cancer and possibly Alzheimer disease.

Antigen: Synthesized peptide derived from the Internal region of human Cathepsin D LC

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in the aorta extracellular space (at protein level).