

**Anti-Human/Mouse/Rat Ephrin A5 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02170**

---

3th Edition

**Description:** May function actively to stimulate axon fasciculation. Induces compartmentalized signaling within a caveolae-like membrane microdomain when bound to the extracellular domain of its cognate receptor. This signaling event requires the activity of the Fyn tyrosine kinase.

**Antigen:** Synthesized peptide derived from the C-terminal region of human Ephrin-A5

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**