

**Anti-Human/Mouse/Rat BCL6 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02230**

---

3th Edition

**Description:** Transcriptional repressor which is required for germinal center formation and antibody affinity maturation. Probably plays an important role in lymphomagenesis.

**Antigen:** Synthesized peptide derived from human Bcl-6 around the non-phosphorylation site of Ser333.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in germinal center T and B cells and in primary immature dendritic cells.