

Anti-Human/Mouse/Rat SMYD2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02264

3th Edition

Description: SET domain-containing proteins, such as SMYD2, catalyze lysine methylation (Brown et al., 2006 [PubMed 16805913]).

Antigen: Synthetic peptide of human SMYD2

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: