

**Anti-Human/Mouse WNT2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02290**

---

3th Edition

**Description:** This gene is a member of the WNT gene family. The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Alternatively spliced transcript variants have been identified for this gene.

**Antigen:** Synthetic peptide of human WNT2

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in brain in the thalamus, in fetal and adult lung and in placenta.