

Anti-Human/Mouse BAZ1B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02315

3th Edition

Description: This gene encodes a member of the bromodomain protein family. The bromodomain is a structural motif characteristic of proteins involved in chromatin-dependent regulation of transcription. This gene is deleted in Williams-Beuren syndrome, a developmental disorder caused by deletion of multiple genes at 7q11.23.

Antigen: Synthetic peptide of human BAZ1B

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitously expressed with high levels of expression in heart, brain, placenta, skeletal muscle and ovary.