

Anti-Human IGFBP6 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02324

3th Edition

Description:IGFBP6 (Insulin Like Growth Factor Binding Protein 6) is a Protein Coding gene. Diseases associated with IGFBP6 include Exfoliative Ichthyosis and Insulin-Like Growth Factor I. Among its related pathways are Metabolism of proteins and Development IGF-1 receptor signaling. GO annotations related to this gene include insulin-like growth factor binding and insulin-like growth factor II binding. An important paralog of this gene is IGFBP5.

Antigen:Synthesized peptide derived from the Internal region of human IGFBP6.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: