

Anti-Human/Mouse C1QL2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02326

3th Edition

Description:C1QL2 (Complement C1q Like 2) is a Protein Coding gene. GO annotations related to this gene include identical protein binding. An important paralog of this gene is C1QL3.

Antigen:Synthesized peptide derived from the C-terminal region of human C1qL2

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: