

**Anti-Human CGRRF1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02340**

---

3th Edition

**Description:**CGRRF1 (Cell Growth Regulator With Ring Finger Domain 1) is a Protein Coding gene.

**Antigen:**Synthesized peptide derived from the C-terminal region of human CGR19.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Ubiquitously expressed with high expression in testis and the cerebellum.