

Anti-Human CGRRF1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02340

3th Edition

Description:CGRRF1 (Cell Growth Regulator With Ring Finger Domain 1) is a Protein Coding gene.

Antigen:Synthesized peptide derived from the C-terminal region of human CGR19.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Ubiquitously expressed with high expression in testis and the cerebellum.