

Anti-Human XKRX Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02346

3th Edition

Description: This gene encodes a protein that is related to a component of the XK/Kell complex of the Kell blood group system. The encoded protein includes several transmembrane domains, is known to be exposed to the cell surface, and may function as a membrane transporter.

Antigen: Synthetic peptide of human XKRX

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed predominantly in the placenta, in syncytiotrophoblasts. Moderate levels in the adrenal gland, low levels in the trachea and very low levels in the bone marrow.