

**Anti-Human SLC10A1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02350**

---

3th Edition

**Description:** The hepatic sodium/bile acid uptake system exhibits broad substrate specificity and transports various non-bile acid organic compounds as well. It is strictly dependent on the extracellular presence of sodium.

**Antigen:** Synthesized peptide derived from the C-terminal region of human NTCP.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**