

Anti-Human SLC10A1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02350

3th Edition

Description:The hepatic sodium/bile acid uptake system exhibits broad substrate specificity and transports various non-bile acid organic compounds as well. It is strictly dependent on the extracellular presence of sodium.

Antigen:Synthesized peptide derived from the C-terminal region of human NTCP.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: