

**Anti-Human CD71 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02356**

---

3th Edition

**Description:**CD71 is a cell surface receptor necessary for cellular iron uptake by the process of receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Mice that are deficient in this receptor show impaired erythroid development and abnormal iron homeostasis.

**Antigen:**Synthesized peptide derived from the N-terminal region of human CD71.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**