

Anti-Human ZNF365 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02376

3th Edition

Description: This gene encodes several isoforms which have different expression patterns and functions. Mutation in this gene is associated with uric acid nephrolithiasis (UAN). Alternatively spliced variants, encoding distinct proteins, have been identified.

Antigen: Synthetic peptide of human ZNF365

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Isoform 4 is expressed in placenta, lung, kidney and pancreas.1 Publication