

Anti-Human/Mouse CARD10 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02436

3th Edition

Description: Activates NF-kappa-B via BCL10 and IKK.

Antigen: Synthesized peptide derived from the Internal region of human CARD 10

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Detected in adult heart, kidney and liver; lower levels in intestine, placenta, muscle and lung. Also found in fetal lung, liver and kidney.