

Anti-Human CD226 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02446

3th Edition

Description: Receptor involved in intercellular adhesion, lymphocyte signaling, cytotoxicity and lymphokine secretion mediated by cytotoxic T-lymphocyte (CTL) and NK cell.

Antigen: Synthesized peptide derived from the Internal region of human DNAM-1.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed by peripheral blood T-lymphocytes.