



**Anti-Human CASP12 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02453**

---

3th Edition

**Description:** Has no protease activity. May reduce cytokine release in response to bacterial lipopolysaccharide during infections. Reduces activation of NF-kappa-B in response to TNF.

**Antigen:** Synthesized peptide derived from the Internal region of human Caspase12.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**