

Anti-Human CASP2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02454

3th Edition

Description:Involved in the activation cascade of caspases responsible for apoptosis execution. Might function by either activating some proteins required for cell death or inactivating proteins necessary for cell survival.

Antigen:Synthesized peptide derived from human Caspase-2 around the non-phosphorylation site of Ser157.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: