

**Anti-Human/Mouse/Rat CASP6 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02456**

---

3th Edition

**Description:**Involved in the activation cascade of caspases responsible for apoptosis execution. Cleaves poly(ADP-ribose) polymerase in vitro, as well as lamins. Overexpression promotes programmed cell death.

**Antigen:**Synthesized peptide derived from human Caspase-6 around the non-phosphorylation site of Ser257.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**