

**Anti-Human/Mouse NCBP2 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02513**

---

3th Edition

**Description:**NCBP2 is a component of the nuclear cap-binding protein complex (CBC), which binds to the monomethylated 5' cap of nascent pre-mRNA in the nucleoplasm. It has an RNP domain commonly found in RNA binding proteins, and contains the cap-binding activity.

**Antigen:**Synthesized peptide derived from the N-terminal region of human CBP20

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**