

**Anti-Human/Mouse/Rat CD239 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02581**

---

3th Edition

**Description:**Laminin alpha-5 receptor. May mediate intracellular signaling.

**Antigen:**Synthesized peptide derived from the Internal region of human CD239.

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Wide tissue distribution (highest in the pancreas and very low in brain). Closely associated with the basal layer of cells in epithelia and the endothelium of blood vessel walls.