

Anti-Human CD223 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02588

3th Edition

Description:Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and chromosomal localization all indicate a close relationship of LAG4 to

Antigen:Synthesized peptide derived from the Internal region of human CD223.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:On cell surface of activated NK and T-lymphocytes.