

**Anti-Human CD70 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02590**

---

3th Edition

**Description:** Cytokine that binds to CD27. Plays a role in T-cell activation. Induces the proliferation of costimulated T-cells and enhances the generation of cytolytic T-cells.

**Antigen:** Synthesized peptide derived from the Internal region of human CD70

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**