



Anti-Human/Mouse/Rat HDGFL3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02594

3th Edition

Description:HDGFL3 (HDGF Like 3) is a Protein Coding gene. An important paralog of this gene is PSIP1.

Antigen:Synthesized peptide derived from the Internal region of human HRP-3.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Detected in testis, heart, spinal cord and brain.