

Anti-Human/Mouse/Rat CDC20 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02612

3th Edition

Description: CDC20 appears to act as a regulatory protein interacting with several other proteins at multiple points in the cell cycle. It is required for two microtubule-dependent processes, nuclear movement prior to anaphase and chromosome separation.

Antigen: Synthesized peptide derived from the Internal region of human Cdc20

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: