

Anti-Human/Mouse/Rat CDK1/2/3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02637

3th Edition

Description: Plays a key role in the control of the eukaryotic cell cycle. It is required in higher cells for entry into S-phase and mitosis. p34 is a component of the kinase complex that phosphorylates the repetitive C-terminus of RNA polymerase II.

Antigen: Synthesized peptide derived from human Cdk1/2/3 around the non-phosphorylation site of Thr14.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: breast cancer tissues, cancer cell lines and glioblastoma tissue.