



Anti-Human FLT3LG Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02642

3th Edition

Description: Stimulates the proliferation of early hematopoietic cells by activating FLT3. Synergizes well with a number of other colony stimulating factors and interleukins.

Antigen: Synthesized peptide derived from the C-terminal region of human Flt3-L.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: