

Anti-Human/Mouse/Rat CDK2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02650

3th Edition

Description: Serine/threonine-protein kinase involved in the control of the cell cycle; essential for meiosis, but dispensable for mitosis.

Antigen: Synthesized peptide derived from human Cdk2/Cdc2 around the non-phosphorylation site of T160.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: