

**Anti-Human/Rat CDKAL1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02656**

---

3th Edition

**Description:**The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.

**Antigen:**Synthesized peptide derived from the N-terminal region of human CDKAL1

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Expressed in pancreatic islets.