

**Anti-Human/Rat ZNF384 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02668**

---

3th Edition

**Description:** Transcription factor that binds the consensus DNA sequence [GC]AAAAA. Seems to bind and regulate the promoters of MMP1, MMP3, MMP7 and COL1A1.

**Antigen:** Synthesized peptide derived from the N-terminal region of human ZNF384.

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**