

Anti-Human/Monkey CEP70 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02681

3th Edition

Description: Plays a role in the organization of both preexisting and nascent microtubules in interphase cells. During mitosis, required for the organization and orientation of the mitotic spindle.

Antigen: Synthesized peptide derived from the Internal region of human CEP70

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: