

Anti-Human/Mouse/Rat FH Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02683

3th Edition

Description:Fumarate hydratase,Also acts as a tumor suppressor

Antigen:Recombinant Protein

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in red blood cells; underexpressed in red blood cells (cytoplasm) of patients with hereditary non-spherocytic hemolytic anemia of unknown etiology.