

**Anti-Human/Mouse CNOT7 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02952**

---

3th Edition

**Description:**The protein encoded by this gene binds to an anti-proliferative protein, B-cell translocation protein 1, which negatively regulates cell proliferation. Binding of the two proteins, which is driven by phosphorylation of the anti-proliferative protein, causes signaling events in cell division that lead to changes in cell proliferation associated with cell-cell contact. The encoded protein downregulates the innate immune response and therefore provides a therapeutic target for enhancing its antimicrobial activity against foreign agents. Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 1 and X.CNOT7 (CCR4-NOT Transcription Complex Subunit 7) is a Protein Coding gene. Among its related pathways are Deadenylation-dependent mRNA decay and Gene Expression. GO annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is CNOT8.

**Antigen:**Synthesized peptide derived from the N-terminal region of human CNOT7

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**