

Anti-Human /Mouse WDR74 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02970

3th Edition

Description: This protein contains six WD repeats. The subcellular location of this protein is Nucleus. Wdr74 is required for RNA transcription, processing and/or stability during preimplantation development and is an essential gene in the mouse.

Antigen: Synthetic peptide of human WDR74

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: